



江苏英科医疗制品有限公司

Jiangsu Intco Medical Products Co.,Ltd

Tel: (0086) 511 83174088 Fax: (0086) 511 83174188

NO.77 Yandunshan Road, Dagang Zhenjiang, Jiangsu Province 212132 P.R. China

Website: www.intcowheelchair.com / www.intcomedical.com / www.intco.com.cn

# SPECIFICATION SHEET



Item No.	YK9132 Steel 18", 20", 22"
Frame	Double cross brace, Powder coating
Front wheel	8*2", PU tire
Rear wheel	24" MAG wheels, PU tire with hand brakes and attendant brakes
Seat width	450mm; 510mm;560mm
Seat depth	400mm
Armrest	Detachable with desk length, PU armpad
Footrest	Elevating footrest
Footplate	Plastic
Seat and backrest	Vinyl upholstery
Weight Capacity (kg)	120kg
G./N. Weight (kg)	29.5/26.6kg
Special	Attendant brakes; Reclining backrest 90°-180° Headrest; Anti-tipper;
Carton size	920*310*960mm
Package	77pcs/20"GP; 224pcs/40"hq

The YK9132 Reclining Wheelchair is intended for patients suffering from low blood sugar or poor circulation. It is made of steel yet can be easily manoeuvred, and is foldable. Front and rear wheels are puncture proof  
Floor to seat Height : 500mm